



## EVA Infusion Bags

- › Made with EVA (ethylene vinyl acetate)
- › 100% PVC-, DEHP-, and latex-free
- › Practical and sturdy carrying handles
- › Each is equipped with three connection points:
  - Luer lock adapter for easy filling and infusing
  - Port for safe medication injection through synthetic septum cap
  - Bag spike adapter for secure connection to infusion set

## Transparent Infusion Bags

- › Available in all standard sizes from 100 to 5000 mL
- › 1000/2000/3000 mL infusion bags with optional three-way transfer set for easy-to-use gravity filling



List Number	Description	Qty
SN2010	EVA Infusion Bag 100 mL	40
SN2020	EVA Infusion Bag 100-300 mL	40
SN2050	EVA Infusion Bag 500 mL	40
SN2100	EVA Infusion Bag 1000 mL	35
SN2200	EVA Infusion Bag 2000 mL	15
SN2300	EVA Infusion Bag 3000 mL	15
SN2400	EVA Infusion Bag 4000 mL	15
SN2500	EVA Infusion Bag 5000 mL	10
SN3100	EVA Infusion Bag 1000 mL three-way transfer set	15
SN3200	EVA Infusion Bag 2000 mL three-way transfer set	10
SN3300	EVA Infusion Bag 3000 mL three-way transfer set	10
SN4020	EVA Infusion Bag 300 mL with 1 access	40
SN4050	EVA Infusion Bag 500 mL with 1 access	40

## Light-Protected Infusions Bags

- › Protection for light-sensitive solutions including lipids, vitamins, and oncology medications
- › Effective UV protection within 89% to 99% of the relevant wavelength range from 175 nm to 525 nm
- › Transparent color allows for visual assessment of the solution
- › Toxicologically harmless material
- › Available in all standard sizes from 100 to 4000 mL
- › Practical and sturdy carrying handles

List Number	Description	Qty
SN2015G	EVA Infusion Bag 150 mL Light-protected	40
SN2020G	EVA Infusion Bag 250 mL Light-protected	40
SN2050G	EVA Infusion Bag 500 mL Light-protected	40
SN2100G	EVA Infusion Bag 1000 mL Light-protected	30
SN2200G	EVA Infusion Bag 2000 mL Light-protected	30
SN2300G	EVA Infusion Bag 3000 mL Light-protected	30



## Multi-Chamber Infusion Bags

- › Secure chamber separation
- › Easy-to-open chamber locks
- › Optional equipment for simultaneous bilateral chamber filling
- › Available three-chamber light-protected bag
- › Protection for light-sensitive solutions including lipids, vitamins, and oncology medications
- › Available in sizes from 500 to 3000 mL
- › Each chamber has a luer lock adapter for easy filling
- › Bags also include a spike connection and a synthetic septum port for injecting medications

List Number	Description	Qty
SN5050	Two-Chamber EVA Infusion Bag 500 mL	30
SN5100	Two-Chamber EVA Infusion Bag 1000 mL	15
SN5200	Two-Chamber EVA Infusion Bag 2000 mL	15
SN5300	Two-Chamber EVA Infusion Bag 3000 mL	15
SN5103G	Three-Chamber EVA Infusion Bag 500/300/200 mL Light-protected	15
SN5203	Three-Chamber EVA Infusion Bag 500/750/750 mL	15



**NEO CARE**  
medical products

Neo Care, a Subsidiary of ICU Medical  
Neo Care GmbH, Altener Str. 136,  
58513 Lüdenscheid  
tel +49 (0) 2351 / 9548-0  
fax +49 (0) 2351 / 9548-20  
email info@neocare.de

**icumedical**  
human connections

800.824.7890 | www.icumed.com